SILI-TEX ORIGINAL

SitOnIt • Seating®

SPECIFICATIONS

PATTERN: Sili-Tex Original

PRICE: FG 4

LEAD TIME: Task/Executive/Side: 2-day Lounge: 10-day

WARRANTY: 5 years

AVAILABILITY: Approved for use on all seating models.

CHARACTERISTICS

Content: 100% silicone Finish: None

Backing: 100% polyester

Width: 54"

Pattern Repeat Size: N/A

Weight Per Linear Yard: 2.2 lbs
Cleaning Code: WS / BC (1:4)
Country of Origin: China

PERFORMANCE

Abrasion: 350,000 D.R.

Seam Slippage (Warp): 182 lbs

Tensile Strength (Warp): 500N +/- 10%

Flammability: CA TB117-2013 Class 1, NFPA

260/UFAC Class 1, ASTM E84

(adhered)

Colorfastness to Light Rating: 4.5
Lightfastness: 200 hours
Crocking, Dry: Class 5
Piling Resistance: N/A

Note: Multiple factors affect fabric

durability and appearance retention, including end-user application and proper maintenance. Wyzenbeek results above 100,000 double rubs have not been shown to be an indicator of increased lifespan.

★ ★ OA

ENVIRONMENTAL

PVC/PU Free HHI Compliant

CAL01350 certified

Qualifies for LEED v4, BREEAM, Well Building, and Living Building Challenge

Advanced solvent free technology

FR Free

California PROP 65 Compliant

No formaldehyde

No plasticizers, phthalates, heavy metals or BPA

No lead Non-toxic

Inherently anti-mildew

No anti-bacterial additives or biocides



SILI-TEX ORIGINAL

SitOnIt • Seating®

CARE + MAINTENANCE GUIDE

STAIN RESISTANCE

Ink: Class 5 Permanent Marker: Class 5 Coffee: Class 5 Blood: Class 5 Class 5 Urine: Ketchup/Mustard: Class 5 Coffee/Wine: Class 5 Class 5 Lipstick: Class 5 Denim:

Rating System

Class 5 - No evidence of stain present

Class 4 – Slight evidence of stain present

Class 3 – Noticeable evidence of stain present

Class 2 - Considerable evidence of stain present

Class 1 – Excessive evidence of stain present

GENERAL MAINTENANCE RECOMMENDATIONS

- Clean and maintain your furniture regularly. Rinse with clean water and dry.
- Clean with a damp cloth with a solution of mild liquid soap and water.
- · A professional cleaning service is recommended.
- · Do not saturate the furniture with liquids.
- Abrasive household cleaners and steel wool should be avoided.
- Treat stains quickly. The stain will become increasingly more difficult to remove the longer it sits.
- Pretest any cleaning method or agent in an inconspicuous area of the furniture before using in a visible location.

STAIN AND DISINFECTIENT RECOMMENDATIONS

For loose dirt and dust:

Clean with a soft sponge or cloth using water.

For stains such as ketchup/mustard or wine/coffee:

Clean with a soft sponge or cloth using water.

For stubborn stains like sunscreen, mascara, soy sauce, pen

ink or permanent marker:

Mix 1 oz of detergent with $\frac{1}{4}$ gallon of water. Use a soft cloth to wash and make sure to rinse the mix off the furniture with water.

For stubborn stains like lipstick or denim dye:

Clean with a soft sponge using 95% ethyl alcohol or 90% Isopropyl alcohol.

For disinfecting:

Use a 1:4 bleach to water mix.

Note:

Make sure to rinse all cleaning agents off the surface with water and a soft cloth. Residue left by cleaning agents will degrade and shorten the life and/or affect the product's appearance if not removed.

APPROVED CLEANERS

This list represents some of the most common options available and is not intended to be comprehensive.

3M HB Quat Disinfectant 25L 3M HB Quat Disinfectant 23H

Accel Prevention Accel TB Wipes

Birex SF

Bleach solution 30%

Cavi Wipes XL

Citrace

Clorox Disinfectant Wipes + Sprays

Clorex Bleach Germicidal Wipes

Clorox Hydrogen Peroxide Cleaner

Dispatch Hospital Cleaner

Diversey Virex II 256 Broad Spectrum Disinfectant

Ecolab 456

Ecolab ASepticare TB+11

Ecolab Discide

Fantastik

Formula 409

Hydrogen Peroxide

Isopropyl Alcohol 90%

Johnson Diversey Disinfectants

Lysol Cleaning Wipes

Lysol Spray

Maxi-Blend

Optim 33 TB Wipes

Oxivir Five 16 Concentrate

Oxivir TB Spray

Oxivir TB Wipes

Oxycide

PDI AF3 Wipes

PDI Plus Wipes

PDI Sani-Cloth

PDI Super Sani-Cloth

Peridox RTU Perisept

Rescue

Soap and Water

Spartan HDQ Neutral

Uniquat Neutral Disinfectant 256

Viragaurd

Virex II 256

Virex TB

Virox 5 Concentrate

Virox 5 RTU

Virox Wipes

SitOnIt • Seating®

